



User's Manual of Safety Goggle EP02

Safety Goggle is non-sterile disposable anti-fog goggle.

1. Brand: Refine
2. Product name: Safety Goggle
3. Material: PVC 70% PC 30%
4. Model: EP02
5. Composition: composed of protective cover and fixing device made of polymer material. Non-sterile, single use.
6. Product Standards: GB 14866-2006 / EN166:2001 / ANSI/ISEA Z87.1-2015
7. Intended use : Safety Goggle is common eye protection equipment, intended to protect against mechanical impact, mechanical injuries, liquids, ocular contact, to prevent crossing infection between medical staffs and patients. Recommended industry: Facility Sanitation, Food Processing, Food Safety, Hospital, Medical Safety.
8. Product instruction: The product and the protective packaging should be checked before using. Stop using it if there is any damage.
9. Precautions:
 - a. One-time use only, No Maintenance, Cleaning required due to being Disposable.
 - b. Use in accordance with the instructions.
 - c. Never use it when damages are found.
 - d. Stay away from chemicals.
 - e. If protection against high speed particles at extremes of temperature is required then the selected eye-protector should be marked with the letter T immediately after the impact letter i.e. FT, BT or AT. If the impact letter is not followed by the letter T then the eye protector shall only be used against high speed particles at room temperature.
 - f. The goggles do not guarantee protection from Harmful Biological Agents.
10. Storage and transport.

















Please store in clean, dry and ventilated indoor place, with relative humidity (10%-93%), Atmospheric pressure for storage (70kPa-106kPa), temperature (-20°C- +40°C), avoid corrosive gas.

Please transport the goggles in original packing.

Warning:

- a. Materials which may come into contact with the wearer's skin could cause allergic reactions to susceptible individuals;
- b. Scratched or damaged oculars should be replaced;
- c. eye-protectors against high speed particles worn over standard ophthalmic spectacles may transmit impacts, thus creating a hazard to the wearer.
- d. These protectors are intended for indoor use where no optical radiation hazards exist.

Symbol instruction

Symbol	Instruction	Symbol	Instruction
	Warning, Caution and Important ! Check the Instruction Manual		Manufacturer
	Date of manufacture		Handle with care
	Recovery		Consult the accompanying documents
	Keep dry		Single-Use
	Atmospheric pressure for storage		Humidity limitation for storage
	Temperature limitation for storage		CE certification
	Lot number		MDEL certification
	EU Representative		According to ANSI/ISEA Z87.1-2015: D3-Special use, used to prevent Droplet / splash
REFINE EP02 1FN-EN166 3 F	<p>According to EN 166:2001, REFINE - manufacturer 1 - Optical class (Class I) F - Protection against high-speed particles, Low energy impact (F) : 45 (+1.5, -0) m/s N - Symbol for Resistance to fogging of oculars EN 166 - Product Standards 3 - Field of use (Liquid Droplet/splash) F - Protection against high-speed particles, Low energy impact (F) : 45 (+1.5, -0) m/s</p>		



Guilin Refine Medical Instrument Co., Ltd.

No.8-3, Information Industrial Park, High-Tech Zone, Qixing District, Guilin, Guangxi, 541004, P.R. China



MedPath GmbH (PPE)
Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany



MedNet EC-REP GmbH (Medical)
Borkstrasse 10, 48163 Münster, Germany

Shelf life: 2 Years